



Eyesight Improvement Device by Color Impulse Therapy

# FORVISION

## CIT-100

Based on color impulse therapy technology and ease of use, the new CIT-100 is ground-breaking device which prevents loss of eyesight and recover eyesight.

By combination of respective illumination level, time, shape, and order regarding frequency and energy of three LED colors (RED, GREEN, BLUE), the light ray, through eye, is absorbed, which makes the best effect for sight improvement. CIT-100 consists of electromagnetic radiation, radiation-illuminating part, light control part, special glasses with emitter of within fields of view.




# FORVISION CIT-100

Product Features

- Color impulse therapy with three LED sources (Red, Green and Blue)
- Eyesight improvement and recovery by color impulse therapy
- Prevention against falling of eyesight
- Beneficial effects of releasing stress, recovering fatigue, deep sleeping, and mental & physical stability

- Slim and ergonomic design for comfortable use
- Portable design for anytime and anywhere use



## How to use


One cycle of 6 weeks			
The first two weeks	The first one week break	The second two weeks	The second one week break
<b>Operating rules</b>			

※ During break week, try to do daily massage of five to ten times around eyes by making circle, which is helpful.

## Eyesight enhancement (Improvement) programs

- 1 Short-sightedness : Morning
- 2 Short-sightedness : Any time
- 3 Short-sightedness : Lunch
- 4 Short-sightedness : Evening
- 5 Astigmatism
- 6 Amblyopia
- 7 Low vision
- 8 Eyestrain
- 9 Spasm of accommodation and false short-sightedness (Pseudomyopia)
- 10 Loss of vision (Presbyopia)


## Examples regarding results of eyesight improvement



### • Source

Professor T.P.T in Russia medicine and science academy and University of Maryland in USA

## Part Name



- **Side Cover** For good airy conditions between internal and external parts
- **Operation Indication LED** LED lights up during operation
- **LCD Window** Display of program number and battery residual quantity
- **Program Button** Program selecting button among built-in programs
- **Power Button** Power on/off button
- **Start Button** Start selected program

## Product Specifications

- **Power supply** - 3.7V, 500mA, Lithium polymer battery
- **Consumption** - 0.6W or less
- **Continuous using time with charged battery** - Approximately 210 minutes
- **Product dimensions** - 190(W) x 100(D) x 70(H)mm
- **Product weight** - Approximately 400g (Except for storage box)
- **Charger** - Rated input AC100 ~ 240V, 50 ~ 60Hz, 0.2A / rated output DC4.2V, 750mA

## Product Package

